CC DISK EASY CoCr

Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Viz návod k použití / Potrebno poštivati uputstva za upotrebu / Consultar las instrucciones de uso / Необходимо учитывать инструкции по применению

Uporabno do / Use by / Gültig bis / Použití / Upotrijebiti do / Usar antes / Использовать до

Kontrolna številka / Batch number / Seriennummer / Seriové číslo / Serijski broj / Numero de serie / Серийный номер

Kataloška številka / Catalogue number / Katalognummer / Katalogové číslo / Kataloški broj / Numero de catalogo / Каталогный номер

Medicinski pripomoček / Medical device / Medizinprodukt / Medicinski proizvod / Zdravotniški pripravek / Producto sanitaria / Медицинское изделие

Edinstveni identifikator pripomočka / Unique device identifier / Einmalige Produktkennung / Jedinствена identifikacija proizvoda / Уникальный идентификатор изделия / Jedinečným identifikátorem prostředku / Identificador único del producto

Proizvajalec / Manufacturer / Hersteller / Výrobce / Proizvodač / Fabricante / Производитель

Datum proizvodnje / Date of manufacture / Herstellungsdatum / Datum výroby / Datum proizvodnje / Fecha de fabricacion / Дата изготовления

SI Navodila za uporabo / GB Directions for use
D Gebrauchsanweisung / ES Instrucciones de uso / CZ Návod k použití
HR Uputa za uporabu / RU Инструкции по использованию

SRB

INTERDENT d.o.o.
Zemunska 22, lok 3
RS-11070 Novi Beograd
T/F: +381/11 217 53 74
www.interdent-bg.com

SK

INTERDENT SK s.r.o.
Za dráhou 21
SK-902 01 Pezínok
T: 0903 418 001
E: interdent@interdent.sk

CZ

INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praha
T: +420/274 783 114
F: +420/274 820 130
E:interdent@interdent.cz

HR

INTERDENT d.o.o.
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385/1 3873 644
F: +385/1 38736 17
E: interdent@interdent.hr

Proizvajalec / Producer / Hersteller
Interdent d.o.o. · SI - 3000 CELJE
Opekarniška cesta 26 · T: +386 (0)3 425-62-00
E: info@interdent.cc · www.interdent.cc

Verzija: 04/2021
Datum: 28.10.2021
Made in Slovenia

SI

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GB

Cobalt-based casted discs, intended for producing prosthetic substitutes in CAD/CAM milling machines for metal-ceramic dental restorations, type 4.

CC DISK EASY CoCr is nickel, beryllium, cadmium and lead free and fulfills the recommendation of the standard EN ISO 22674 and EN ISO 9693 for dental alloys.

The alloy has excellent corrosion resistance and good milling properties. It is easy to drill in and gives perfectly polished surface. Only few oxides are formed during ceramic firing.

CC DISK EASY CoCr is a medical device according to European directive 93/42/EEC and is registered by the manufacturer under CE 0197.

Composition (%):	Technical data
Co 62,5	Type 4
Cr 27,2	Temperature solidus, liquidus 1380 °C, 1450 °C
W 8,2	Vickers hardness HV 10 249
Si 1,7	Coefficient of thermal expansion 25 - 500 °C 14,4 x 10 ⁻⁶ K ⁻¹
Mn <1,0	20 - 600 °C 14,6 x 10 ⁻⁶ K ⁻¹
	Elongation limit Rp 0,2 380 MPa
	Tensile strength Rm 553 MPa
	Modul of elasticity (GPa) 167 GPa
	Ductile yield A5 16,9 %

Ceramic firing: You can use all types of ceramic which you also use in conventional method on cast crowns and bridges (VITA, Ivoclar, Ducera, etc.).

Please be careful about cooling procedure. Cobalt-based alloys must be cooled down slowly to 600 °C after each firing step (dentin, correction, glaze).

- Remove connectors and smoothen the sharp edges. Sand-blast the surface with 250 microns aluminum oxide on 2-3 bar air pressure (REF 404 Interlox 250 micr).
- Clean the surface in distilled water with ultrasound, degrease it with ethyl acetate or under vapor.
- Carry out oxide firing (only to control the surface) 5 min at 960°C. After oxide firing removing of oxides is required - use aluminum oxide 250 microns and degrease the surface again.

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- For a good bond between metal and ceramic it is mandatory to use a "metal-ceramic bond". VITA NP BOND paste should be applied in a thin layer, with brush (as Wash opaque). For firing in the ceramic furnace check the below program.**

Firing program using Vita NP Bond:						
Vt. in °C	→ min	↑ min	↓ °C/min	Temp. cca °C	→ min	VAC min
600	6.00	6.00	60	960	1.00	6.00

- When "metal-ceramic bond" phase is finished, proceed with opaque layering, to cover the metal shade of the construction.
- For firing procedure, consult ceramic manufacturer's instructions for use.**
- After each firing step (dentine bake, correction bake, glaze firing), cool down slowly to 600 °C.

Finishing: After last firing step finish the visible surface of the alloy with polishing rubbers and polishing paste (REF 461 Diamond polishing paste) to high gloss. On the cervical parts of the crowns, remove the oxides with glass pearls or very carefully with 50 microns aluminum oxide (REF 406 Interlox 50 micr.) Be careful not to damage the gingival fit!

Side effects: Consider allergic hypersensitivities to contents of the alloy.

Safety precaution! Wear protective clothing, respiratory mask with filter FFP2, gloves, goggles and use suction units when finishing and sand-blasting the construction.

Reciprocal Actions! In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Warranty: Our recommendations for use are based upon our own experience and trials and can only be considered a standard value. Our products are subject to further development. Therefore, alterations in construction and composition are reserved.

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