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Cobalt-based dental casting alloy for partial dentures without beryllium, nickel, cadmium and lead, type 5.

I-MG HE is biocompatible alloy, with a high degree of elasticity for repeated activation of the clasps. It has outstanding resistance to corrosion and high strength. **I-MG HE** can be well polished and laser welded. It meets the requirements of standard EN ISO 22674.

| Composition | (Mass-%) | Properties | |
|-------------|----------|-------------------------------|--|
| Co | 63,5 | Type | 5 |
| Cr | 29,1 | Density | 8,2 g/cm ³ |
| Mo | 5,0 | Solidus, liquidus temperature | 1300 °C, 1330 °C |
| Si | 1,3 | Casting temperature | 1450 °C |
| Mn, C, N | < 1% | Vickers hardness | HV 10 cca. 355 |
| | | Tensile strength | Rm 941 MPa (N/mm ²) |
| | | Yield stress | Rp 0,2 687 MPa (N/mm ²) |
| | | Modulus of elasticity | E cca. 200.000 MPa |
| | | Elongation | A5 6,9 % |

Indications: Used for frames in removable prosthetic dentures, for appliances in which parts require the combination of high stiffness and proof stress, e.g. thin removable partial dentures, parts with thin cross-sections, clasps, attachments, bars.

Contraindications: in case of known allergies to alloy content

Patient selection criteria: Totally or partially edentulous patients.

Intended users: Product is intended to be used by professionals – dental technicians.

Recommendation for use:

Wax up: **Sprues** should be with round profile **Ø 3.5 – 4 mm**, avoid rectangular placing of sprues and direct slap of alloy.

Investing: Use phosphate-bonded partial-denture investment material – **Modelcast** or **Modelcast S**. Preheating temperature of investment is **950 – 1000 °C** and depend on moulding and construction of sprues. At total plates is the end temperature 1050 °C. Holding time of end temperature is **45 – 60 minutes**, and depend on the size of the casting flask and quantity of flasks in the furnace. We ask you to follow recommendation for use of casting machine manufacturer for casting. **Casting:** For

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I-MG HE use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use. When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around ingots. Start casting as soon as the bath begins to vibrate. Do not overheat the alloy. After cooling down to ambient temperature deflask the cast and sand-blasted it under 4 bar pressure with **Interlox 250 µm**. Be careful with clasp sides and stress breakers.

Polishing: For finishing use carbide, ceramically bonded stones or hart metal milling tools.

Use polishing liquid for electrolyte polishing in dental laboratory. Clasps and fitting parts should be protected before polishing with special varnish that prevent uncontrolled material removing.

After treating with stones, milling tools and measuring the fit, the cast should be polished with rubber polisher, and at the end with **Universal polishing paste for Cr-Co-Mo alloys REF 460** up to high gloss.

Soldering and welding: Soldering with Co based lot and high temperature flux. Laser welding with suitable base-metal welding wires. The requirements of EN ISO 9333 must be followed.

Safety warnings: Use suction unit, wear gloves, googles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

Reciprocal Actions! In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Side effects: Consider allergic hypersensitivities to contents of the alloy.

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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