SL ŽICA ZA LASERSKO VARJENJE

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GB

LASER WELDING WIRE

Contains no nickel, beryllium, cadmium and lead and fulfills the recommendations of the standards EN ISO 28319.

Mass composition in %: Co 65,5, Cr 28,0, Mo 5,5, Si, Mn, C

Welding:

Sand blast or grind welding joint in order to remove all oxides or residuals and prepare proper welding gap with tapered walls. Fix the parts i.e. on a cast die, and fill in the Laser Welding Wire beginning on the ground in several layers one onto each other to form a proper bevel weld seam. Do not weld seam at once but alternatively and crosswise from all sides in order to minimize distortion. For proper handling of the wire position just the tip of the wire into laser beam slightly above the surface of parts to be jointed (Material should drop down from tip of the wire to welding zone point by point.) In order to avoid cracks in the welding zone make sure to finish the welding procedure in short time not to have the already welded zone repeatedly heated too often.

Finishing:

Sand-blast the visible surface of the frame with Interlox 50 µm and then polish it with suitable grinding and polishing instruments for nickel alloys. At the end use Universal polishing paste for Co-Cr-Mo alloys REF 460 and polish up to high gloss. Clean polished surface with vapour or boiling water.

Final Safety Notes!

Metal dusts are harmful. Use dust extractor, wear gloves, goggles and protective mask with filter FFP2 during handling with alloy.

Side effects:

Consider allergic hypersensitivities to contents of the alloy.

Reciprocal Actions!

In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Warranty!

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Consultar las instrucciones de uso / Необходимо учитывать инструкции по применению

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