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## WORKING INSTRUCTIONS MODELCAST

Modelcast is temperature resistant phosphate investment material used for making partial dentures with metal base for conventional (by step) heating method. It can be duplicated in silicone or in gel.

- For good mixture, good chemical binding and smooth cast the investment should be mixed 15 seconds by hand and 60 seconds by vacuum.
- Before you mix the powder and liquid the mixing bowl has to be clean and dry.
- When you use the silicon for duplicating and you place it under the pressure, you must continue the complete process under pressure.
- Use different mixing bowl for investment material and plastery. Take care that the investment does not come into the contact with plaster.
- We recommend using degreasing agent for wax parts. (Interwaxit). Interwaxit is necessary well-dry after use.
- Keep furnace door closed for 10 minutes after insertion of the mould! Danger of explosion!

**WORKING INSTRUCTIONS:**

- First pour liquid and then add powder.
- Manual mixing time 15 s
- Vacuum mixing time 60 s
- Duplicate model in silicone or gel.
- Prepare the duplicated model (The base of the model should have a minimum thickness of 1 cm)
- Remove the model from silicon after 30-40 minutes.
- Dip/Spray the SILICONE duplicated model: MODELHART SPRAY
- Dip/Spray the GEL duplicated model: MODELHART, EKOHART
- Do the wax modelation
- Invest into the mould (see the mixing ratio)
- Scratch the surface of the mould before place it in the furnace

**MIXING RATIO POWDER/ LIQUID (DUPLICATION IN GEL)**

200 g / 30 ml.....model

400 g / 60 ml.....mould

	Exposol	Distilled water
Partial denture	Model: 60% (18 ml)	Model: 40% (12 ml)
Combi technique	Mould: 60% (36 ml)	Mould: 40% (24 ml)

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**MIXING RATIO POWDER/ LIQUID (DUPLICATION IN SILICONE)**

200 g / 31ml.....model

400 g / 62 ml.....mould

	Exposol	Distilled water
Partial denture	Model: 65% (20 ml) <p>Mould: 65% (40 ml)</p>	Model: 35% (11 ml) <p>Mould: 35% (22 ml)</p>
Combi technique	Model: 75% (23 ml) <p>Mould: 75% (46 ml)</p>	Model: 25% (8 ml) <p>Mould: 25% (16 ml)</p>

**CONVENTIONAL HEATING METHOD:**

Temperature	Rate of heating	Holding time
23 – 270°С	5°С/min	40 min
270°С – 580°С	7°С/min	30 min
580°С – 950°С	9°С/min	60 min

Holding time minimum .60 minutes – (depending on size and number of moulds).

Casting (I-MG, I-MG FH, I-MG EKO)

Let the mould cool down to the room temperature, do not force the cooling.

After cooling to the room temperature take the cast out of the investment

**WARRANTY!**

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experi-ence and trials and can only be considered as standard values.

Our products are subject to further development. Therefore changes in construction and composition are reserved.

**STORAGE:** at room temperature

**INVESTING TEMPERATURE:** 19°С – 23°С

**PERSONAL PROTECTION:** We recommend wearing half facemask with filter FFP1.

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